

InSpec™

Guide to surface contamination control

Product description	Manufacturer code	Art. No.	Colour code	Chemicals	Mixture	Filtered (0.2 µm)	Bag in bottle	Screw cap	Wipe/ Mop cover	Sterile	Gamma-irradiated in case of wipes	Aseptically treated	Shelf life from date of manufacture	Recommended shelf life after opening		
InSpec™ IPA 500ml Sterile Trigger Spray	IPAWFI-500MLS	8004013	Blue	Isopropanol (Propan-2-ol)	70 30 v/v 70% isopropanol mixed with 30% WFI	✓	✓			✓			2 years	3 months		
InSpec™ IPA 1L Sterile Trigger Spray	IPAWFI1-1LS	8004012				✓	✓					✓			2 years	3 months
InSpec™ IPA 5L Screw cap	IPAWFI1-5LS	8004014				✓		✓				✓			2 years	10 hours
InSpec™ IPA 15 Pouch Wipes	IPAWFIPCH15S	8005010				✓			✓			✓			2 years	10 hours
InSpec™ IPA 50 Pouch Wipes	IPAWFIPCH50S	8005011				✓				✓		✓			2 years	10 hours
InSpec™ IPA Sterile Tub Wipes – 130 wipes / roll	IPAWFI12TB130S	8005012				✓					✓	✓			2 years	n. s.
InSpec™ IPA Tub Wipes – 200 wipes / roll	IPAWFI12TB200S	8005013				✓						✓	✓		2 years	n. s.
InSpec™ IPA Polyester 15 Wipes	IPAWFI15POL9	8005019				✓						✓	✓		2 years	10 hours
InSpec™ IPA Polyester 20 Wipes	IPAWFI20POL12	8005027				✓						✓	✓		2 years	10 hours
InSpec™ IPA Polyester 30 Wipes	IPAWFI30POL12	8005028				✓						✓	✓		2 years	10 hours
InSpec™ IPA Sachet – Single Wipes	IPAWFICH1S	8005014				✓						✓	✓		2 years	2 hours
InSpec™ DE 500ML Sterile Trigger Spray	DEWFI11-500MLS	8004032				Dark Blue	Ethanol, denatured by addition of isopropanol	70 30 v/v 70% andhanol mixed with 30% WFI	✓	✓			✓			2 years
InSpec™ DE 1L Sterile Trigger Spray	DEWFI10-1LS	8004030	✓	✓								✓			2 years	3 months
InSpec™ DE 5L	DEWFI10-5LS	8004031	✓		✓							✓			2 years	10 hours
InSpec™ DE 15 Pouch Wipes	DEWFIPCH15S	8005030	✓							✓		✓			2 years	10 hours
InSpec™ DE 50 Pouch Wipes	DEWFIPCH50S	8005031	✓								✓	✓			2 years	10 hours
InSpec™ DE Sterile Tub Wipes – 130 wipes / roll	DEWFI12TB130S	8005032	✓									✓	✓		2 years	n. s.
InSpec™ DE Sterile Tub Wipes – 200 wipes / roll	DEWFI12TB200S	8005033	✓									✓	✓		2 years	n. s.
InSpec™ OX 1L Sterile Trigger Spray	OXPW40-1LS	8004020	Red	Sporicidal mixture of hydrogen peroxide and peracetic acid	6% hydrogen peroxid / 0.4% peracetic acid mixed with PW				✓	✓			✓		✓	2 years
InSpec™ OX 5L Screw Cap	OXPW40-5LS	8004021				✓		✓				✓		✓	2 years	10 hours
InSpec™ OX Burstable Pouch – 15 Wipes	OXPW40BP15S	8005020				✓				✓		✓		✓	2 years	2 hours
InSpec™ OX Burstable Mop – 2 Sterile Mops	OXPWMOP2S	8005021				✓					✓	✓		✓	2 years	2 hours
InSpec™ QT 900MLS Sterile Trigger Spray	QTWFI30-1LS	8004051	Green	Quaternary ammonium compound	4% quaternary ammonium compound mixed with WFI	✓	✓			✓	✓		2 years	3 months		
InSpec™ QT 5L	QTWFI30-5LS	8004062				✓		✓				✓	✓		2 years	10 hours
InSpec™ QT Concentrate	QTCNC30-100	8004050				✓				✓		✓	✓		2 years	10 hours
InSpec™ QT Burstable Pouch – 15 Wipes	QTWFI30BP15S	8005050				✓					✓	✓	✓		2 years	2 hours
InSpec™ QT Burstable Mop – 2 Sterile Mops	QTFWIMOP2S	8005051				✓					✓	✓	✓		2 years	2 hours
InSpec™ QT+ Sterile 1L Trigger RTU	QTPLWFI31-900MLS	8004061				✓	Quaternary ammonium compound with complexing agents	2% quaternary ammonium compound mixed with WFI	✓	✓			✓	✓		2 years
InSpec™ QT+ Sterile 5L RTU	QTPLWFI31-5LS	8004062		✓		✓						✓	✓		2 years	10 hours
InSpec™ QT+ Sterile Concentrate 102ml	QTPLCNC-102	8004060		✓						✓		✓	✓		2 years	10 hours
InSpec™ QT+ Burstable Wipes – 15 Wipes	QTPLWFI31BP15S	8005060		✓							✓	✓	✓		2 years	2 hours
InSpec™ AN 900MLS Sterile	ANWFI20-1LS	8004071		Orange	Diamine with water	4% diamine mixed with WFI	✓	✓			✓	✓		2 years	3 months	
InSpec™ AN 5L Screw Cap	ANWFI20-5LS	8004072	✓					✓				✓	✓		2 years	10 hours
InSpec™ AN Concentrate	ANCNC20-100	8004070	✓							✓		✓	✓		2 years	10 hours
InSpec™ AN Burstable Pouch – 15 Wipes	ANWFI20BP15S	8005070	✓								✓	✓	✓		2 years	2 hours
InSpec™ AN Burstable Mop – 2 Sterile Mops	ANWFIMOP2S	8005071	✓								✓	✓	✓		2 years	2 hours
InSpec™ HA 1L Sterile Trigger Spray	HAPW50-1LS	8004040	Purple				Sporicidal mixture made of hypochlorous acid	7% hypochlorous acid mixed with PW	✓	✓			✓		✓	1 year
InSpec™ HA 5L	HAPW50-5LS	8004041		✓		✓						✓	✓		1 year	10 hours
InSpec™ HA Burstable Wipes – 15 Wipes	HAPW50BP15S	8005040		✓							✓	✓	✓		1 year	2 hours
InSpec™ HA Burstable Mop – 2 Sterile Mops	HAPWMOP2S	8005041		✓							✓	✓	✓		1 year	2 hours

n. s. = not specified
 RTU = ready to use
 WFI = water for injection
 PW = purified water

Corresponding no-sterile versions available!

Manufacturer's data